

# SCOUTTUBE

## THE HYGIENIC MOUTHPIECE FOR SINGLE-PATIENT USE

ScoutTubes ensure a safe, accurate and patient-friendly spirometry test with the GANSHORN SpiroScout and the SpiroScout SP plus devices.

### The main advantages of the ScoutTubes are:

- ✦ Always accurate
- ✦ Resistance-free
- ✦ Calibration-free
- ✦ Disposable, therefore no cleaning
- ✦ Maintenance-free
- ✦ No errors caused by miscalibration

- ✦ Ensures a tight seal around the lips
- ✦ Designed for all mouth shapes and sizes
- ✦ Suits all ages from children to adults, at all levels of fitness



### POSITIONING OF TEETH

- ✦ Keeps ScoutTube securely in place

### TIGHT LIP SEAL


- ✦ Prevents leaks

### MOUTHPIECE

- ✦ Oval in shape for optimal lip seal and teeth positioning
- ✦ Conical to suit different mouth sizes

### FILTER MEMBRANE

- ✦ Bacteria filter membrane built into each ScoutTube
- ✦ Prevents cross-contamination





# Benefits of ScoutTube

## ECONOMICAL

- ∴ No running costs
- ∴ No costs for calibration equipment
- ∴ No costs or downtime for maintenance

## SIMPLE APPLICATION

- ∴ Just insert the ScoutTube, ready for measurement

## MAINTENANCE-FREE

- ∴ No maintenance needed with ScoutTube or SpiroScout
- ∴ No calibration or checking procedures required

## HYGIENIC

- ∴ One of the safest respiratory flow measurement systems available
- ∴ Can be fitted with the GANSHORN PFT filter

## ENVIRONMENTALLY SAFE

- ∴ Lightweight PE, environmentally friendly, disposable, and bio-degradable



## DISPOSABLE AND SUSTAINABLE

To ensure a sustainable product, SCHILLER compensates the CO<sub>2</sub> emitted during production, packaging, shipment, and disposal of the ScoutTubes and the PFT filters.



ScoutTube is compatible with the following devices:

- ∴ GANSHORN SpiroScout
- ∴ SpiroScout SP plus
- ∴ SpiroScout SP (discontinued)



**SCHILLER**  
The Art of Diagnostics

Headquarters: SCHILLER AG, Altgasse 68, CH-6341 Baar, Switzerland | Phone +41 41 766 42 42 | sales@schiller.ch | www.schiller.ch

SCHILLER reserves the right to make changes in design, specifications and features shown or discontinue the product described at any time without notice or obligation. Images may be representative. Some products may not be regulatory cleared or released for sale in all markets. Contact your SCHILLER representative for the most current information and availability.

Article no.: 2.500743  2023/10, EN | © SCHILLER. ALL RIGHTS RESERVED.

